

■ Made
■ in
■ Germany



The **OTS 40** is a stable and robust slide stretching table that is available in five different sizes and finds its place in every laboratory.

The work surface consists of a plastic-coated aluminum plate that is 6 mm thick. This enables rapid heating to the desired working temperature, which is then kept stable during the entire duration of use. Temperatures of max. +99 °C can be set using a digital temperature controller at the front. The current temperature is shown on the display during operation.

The **OTS 40** series has proven itself particularly thanks to its high temperature constancy and long life time in routine laboratory work.



Technical Data OTS 40.1520		Cat. No.
Temperature:	+30 °C to +99 °C	
Dimensions Overall (W/D/H):	150 x 200 x 80 mm	
Weight:	1.7 kg	
Power Supply:	230 V / 50 – 60 Hz / 120 VA	01-4001-00
	115 V / 50 – 60 Hz / 120 VA	01-4101-00
Technical Data OTS 40.1540		Cat. No.
Temperature:	+30 °C to +99 °C	
Dimensions Overall (W/D/H):	150 x 400 x 80 mm	
Weight:	2.9 kg	
Power Supply:	230 V / 50 – 60 Hz / 300 VA	01-4002-00
	115 V / 50 – 60 Hz / 300 VA	01-4102-00
Technical Data OTS 40.2025		Cat. No.
Temperature:	+30 °C to +99 °C	
Dimensions Overall (W/D/H):	200 x 250 x 80 mm	
Weight:	2.4 kg	
Power Supply:	230 V / 50 – 60 Hz / 250 VA	01-4003-00
	115 V / 50 – 60 Hz / 250 VA	01-4103-00
Technical Data OTS 40.2530		Cat. No.
Temperature:	+30 °C to +99 °C	
Dimensions Overall (W/D/H):	250 x 300 x 80 mm	
Weight:	3.6 kg	
Power Supply:	230 V / 50 – 60 Hz / 350 VA	01-4004-00
	115 V / 50 – 60 Hz / 350 VA	01-4104-00
Technical Data OTS 40.3040		Cat. No.
Temperature:	+30 °C to +99 °C	
Dimensions Overall (W/D/H):	300 x 400 x 80 mm	
Weight:	5.4 kg	
Power Supply:	230 V / 50 – 60 Hz / 550 VA	01-4005-00
	115 V / 50 – 60 Hz / 550 VA	01-4105-00



Distributed by:
Axlab A/S
 +45 35 43 18 81
 axlab@axlab.dk
 www.axlab.dk

Axlab

YOUR VISION.
 OUR INSPIRATION

MEDITE
MEDITE Medical GmbH

Wollenweberstr. 12 • 31303 Burgdorf • Germany

Phone +49 5136 88 84-0 • Fax +49 5136 88 84-55

info@medite.de • www.medite.de

Subject to change without prior notice. Pictures and colours may vary.