

■ Made  
■ in  
■ Germany



The round Tissue Floatation Bath **TFB 45** is a compact and inexpensive device made of powder-coated aluminum with a capacity of 1.5 liters. Paraffin sections can be stretched and drawn onto the slide. The deep black plastic coating provides optimal contrast to the paraffin section and guarantees a long shelf life.

Thanks to its low weight and compact size, it can be used flexibly at various workplaces in pathological, histological and other laboratories.

The Tissue Floatation Bath has a heated rim that is suitable for depositing and gently drying the slides that have been drawn up. For more storage space on a particularly wider rim, the **TFB 55** is available.

With the user-friendly controller with temperature scale, the desired value can be set exactly to  $\pm 2^\circ \text{C}$ . The rim is heated parallelly. Overheating protection is integrated.

Suitable thermometers are available as accessories, as well as an aluminum lid that ensures temperature stability and protects the bath against dirt when it is not in use.



Technical Data TFB 45		Cat. No.
Bath Capacity:	max. 1.5 liters	
Temperature:	+30 °C to +90 °C	
Rim Width:	45 mm	
Dimensions Bath Internal (Ø/H):	210 x 55 mm	
Dimensions Overall (Ø/H):	300 x 90 mm	
Weight:	ca. 2 kg	
Power Supply:	230 V / 50 – 60 Hz / 250 VA	01-4500-00
	115 V / 50 – 60 Hz / 250 VA	01-4510-00
Technical Data TFB 55		
Bath Capacity:	max. 1.5 liters	
Temperature:	+30 °C to +90 °C	
Rim Width:	55 mm	
Dimensions Bath Internal (Ø/H):	210 x 55 mm	
Dimensions Overall (Ø/H):	340 x 90 mm	
Weight:	ca. 2 kg	
Power Supply:	230 V / 50 – 60 Hz / 250 VA	01-5500-00
	115 V / 50 – 60 Hz / 250 VA	01-5510-00


Accessories	Purchasing Unit	Cat. No.
Accessories Kit for <b>TFB 45</b> incl. lid, thermometer and matching holder	1 kit	11-4500-00
Accessories Kit for <b>TFB 55</b> incl. lid, thermometer and matching holder	1 kit	11-5500-00
Thermometer for <b>TFB 45</b> and <b>TFB 55</b> , 0 – 100 °C	1 pc.	51-4551-00
Lid made of aluminum for <b>TFB 45</b> and <b>TFB 55</b>	1 pc.	11-4540-00
Thermometer Holder for <b>TFB 45</b>	1 pc.	11-4541-00
Thermometer Holder for <b>TFB 55</b>	1 pc.	11-5541-00



Distributed by:  
**Axlab A/S**  
 +45 35 43 18 81  
 axlab@axlab.dk  
 www.axlab.dk

**Axlab**

YOUR VISION.  
 OUR INSPIRATION



**MEDITE Medical GmbH**

Wollenweberstr. 12 • 31303 Burgdorf • Germany  
 Phone +49 5136 88 84-0 • Fax +49 5136 88 84-55  
 info@medite.de • www.medite.de

Subject to change without prior notice. Pictures and colours may vary.

